



## Enteral Safety System *by NeoMed, Inc.*


### NEOMED PROUDLY ANNOUNCES THE FIRST COMPREHENSIVE ENTERAL SAFETY SYSTEM (ESS) AVAILABLE ON THE MARKET TODAY.

Our Enteral Safety System includes oral/enteral syringes, enteral only extension sets and feeding tubes that are designed to be incompatible with any luer or IV devices. Additionally, NeoMed's Enteral Safety System (ESS) complies with all recommendations set forth by The Joint Commission Journal on Quality and Patient Safety, May 2008, Vol. 34 Num. 5 pg. 285-292, "Enteral Feeding Misconnections: A Consortium Position Statement" and The Journal of Parenteral and Enteral Nutrition Vol. 20 Num. 10 pg. 1-46, "Enteral Nutrition Practice Recommendations." By utilizing NeoMed's ESS, the risk of non-compliance and the occurrence of unintended consequences is functionally ELIMINATED.

At the cornerstone of the Enteral Safety System (ESS) is our full line of oral/enteral dispensers (1ml, 3ml, 6ml, 12ml, 20ml, 35ml and 60ml) that facilitate pump feedings and nutrition delivery through our enteral only extension sets down to our enteral only feeding tubes. For your convenience, we listed the product item numbers, descriptions and packaging information for all of our dispensers, extension sets and feeding tubes.

NeoMed offers the only complete enteral delivery system that complies with all recommendations set by "Dangerous Links Persist Between Incompatible Tubing System." Inside the Joint Commission Online, Vol. 13 No.12 pg. 1-8.

**Please see back for ordering information**



Orange graduation markings and text on the barrel of the syringe are a ®registered trademark of NeoMed, Inc.

## Product List

Product Code	Product Description	Quantity
<b>NeoMed Silicone Feeding Tubes</b>		
FTS5.0S-EO	Enteral F/T Radiopaque Orange Stripe 5FrX40cm	10/Box
FTS6.5S-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX40cm	10/Box
FTS8.0S-EO	Enteral F/T Radiopaque Orange Stripe 8FrX40cm	10/Box
FTM5.0S-EO	Enteral F/T Radiopaque Orange Stripe 5FrX60cm	10/Box
FTM6.5S-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX60cm	10/Box
FTM8.0S-EO	Enteral F/T Radiopaque Orange Stripe 8FrX60cm	10/Box
FTL5.0S-EO	Enteral F/T Radiopaque Orange Stripe 5FrX90cm	10/Box
FTL6.5S-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX90cm	10/Box
FTL8.0S-EO	Enteral F/T Radiopaque Orange Stripe 8FrX90cm	10/Box
<b>NeoMed Polyurethane Feeding Tubes</b>		
FTS5.0P-EO	Enteral F/T Radiopaque Orange Stripe 5FrX40cm	10/Box
FTS6.5P-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX40cm	10/Box
FTS8.0P-EO	Enteral F/T Radiopaque Orange Stripe 8FrX40cm	10/Box
FTM5.0P-EO	Enteral F/T Radiopaque Orange Stripe 5FrX60cm	10/Box
FTM6.5P-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX60cm	10/Box
FTM8.0P-EO	Enteral F/T Radiopaque Orange Stripe 8FrX60cm	10/Box
FTL5.0P-EO	Enteral F/T Radiopaque Orange Stripe 5FrX90cm	10/Box
FTL6.5P-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX90cm	10/Box
FTL8.0P-EO	Enteral F/T Radiopaque Orange Stripe 8FrX90cm	10/Box
<b>Orange Enteral Extension Sets</b>		
NM-60ENEN	60" Orange male/female Enteral Only set	50/case
IC-40ENENB	40" Orange male/female Enteral Only set with extra large bore	50/case
IC-60ENENY	60" Orange male/female Enteral Only set with Enteral Only Y-site	50/case
<b>Individually Packaged Dispensers (Pump Compatible)</b>		
NM-S1EO	Individually Packaged Oral / Enteral Dispensers 1ml	200/Box
NM-S3EO	Individually Packaged Oral / Enteral Dispensers 3ml	200/Box
NM-S6EO	Individually Packaged Oral / Enteral Dispensers 6ml	200/Box
NM-S12EO	Individually Packaged Oral / Enteral Dispensers 12ml	200/Box
NM-S20EO	Individually Packaged Oral / Enteral Dispensers 20ml	100/Box
NM-S35EO	Individually Packaged Oral / Enteral Dispensers 35ml	100/Box
NM-S60EO	Individually Packaged Oral / Enteral Dispensers 60ml	100/Box
<b>Sterile Syringes (Pump Compatible)</b>		
SD-S1EO	Sterile Oral / Enteral Syringes 1ml	200/bx
SD-S3EO	Sterile Oral / Enteral Syringes 3ml	200/bx
SD-S6EO	Sterile Oral / Enteral Syringes 6ml	200/bx
SD-S12EO	Sterile Oral / Enteral Syringes 12ml	200/bx
SD-S20EO	Sterile Oral / Enteral Syringes 20ml	100/bx
SD-S35EO	Sterile Oral / Enteral Syringes 35ml	100/bx
SD-S60EO	Sterile Oral / Enteral Syringes 60ml	100/bx



Orange signals Enteral Safety.  
Orange on the barrels signals authentic NeoMed.

NeoMed also carries a full line of pharmacy dispensers.

4-0810