



Polyurethane Feeding Tubes *by NeoMed, Inc.*


SINCE NEOMED'S INCEPTION, DELIVERING ENTERAL FEEDINGS SAFELY HAS BEEN THE HALLMARK OF OUR PRODUCTS.

Until now, feeding tube design has largely focused on creating hub connections that avoid dangerous misconnection errors. The NeoMed Polyurethane Feeding Tube line moves beyond hub design to address factors that affect patient safety and outcome.

It is our philosophy that feeding tubes should be designed to reduce tissue trauma common with oversized and sharp side port holes, or sharp unfinished distal tips. Bacteria growth can be attributed to dangerous material reservoirs within closed distal tips or bifurcated hub designs.

NeoMed's Polyurethane Feeding Tube features and benefits include:

- Open distal tip with no sharp edges or hidden cavities, eliminates residual colonization associated with closed tip feeding tubes
- Insertion markings at each centimeter on the radiopaque tubing help to confirm placement
- Soft medical grade polymer supports quality developmental care
- Indwell time up to 30 days minimizing additional tube placements
- Provides a cost effective solution



NeoMed also offers a full line of Oral Syringes

Please see back for ordering information

Product List

Product Code	Product Description	Quantity
NeoMed Polyurethane Feeding Tubes (Enteral Only Connector)		
FTS5.0P-EO	Enteral F/T Radiopaque Orange Stripe 5FrX40cm	10/Box
FTS6.5P-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX40cm	10/Box
FTS8.0P-EO	Enteral F/T Radiopaque Orange Stripe 8FrX40cm	10/Box
FTM5.0P-EO	Enteral F/T Radiopaque Orange Stripe 5FrX60cm	10/Box
FTM6.5P-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX60cm	10/Box
FTM8.0P-EO	Enteral F/T Radiopaque Orange Stripe 8FrX60cm	10/Box
FTL5.0P-EO	Enteral F/T Radiopaque Orange Stripe 5FrX90cm	10/Box
FTL6.5P-EO	Enteral F/T Radiopaque Orange Stripe 6.5FrX90cm	10/Box
FTL8.0P-EO	Enteral F/T Radiopaque Orange Stripe 8FrX90cm	10/Box
Orange Enteral Extension Sets		
NM-60ENEN	60" Orange male/female Enteral Only set	50/Case
IC-40ENENB	40" Orange male/female Enteral Only set with extra large bore	50/Case
IC-60ENENY	60" Orange male/female Enteral Only with Enteral Only Y-Site	50/Case
Individually Packaged Dispensers (Pump Compatible)		
NM-S1EO	Individually Packaged Oral / Enteral Dispensers 1ml	200/Box
NM-S3EO	Individually Packaged Oral / Enteral Dispensers 3ml	200/Box
NM-S6EO	Individually Packaged Oral / Enteral Dispensers 6ml	200/Box
NM-S12EO	Individually Packaged Oral / Enteral Dispensers 12ml	200/Box
NM-S20EO	Individually Packaged Oral / Enteral Dispensers 20ml	100/Box
NM-S35EO	Individually Packaged Oral / Enteral Dispensers 35ml	100/Box
NM-S60EO	Individually Packaged Oral / Enteral Dispensers 60ml	100/Box
Sterile Syringes (Pump Compatible)		
SD-S1EO	Sterile Oral / Enteral Syringes 1ml	200/bx
SD-S3EO	Sterile Oral / Enteral Syringes 3ml	200/bx
SD-S6EO	Sterile Oral / Enteral Syringes 6ml	200/bx
SD-S12EO	Sterile Oral / Enteral Syringes 12ml	200/bx
SD-S20EO	Sterile Oral / Enteral Syringes 20ml	100/bx
SD-S35EO	Sterile Oral / Enteral Syringes 35ml	100/bx
SD-S60EO	Sterile Oral / Enteral Syringes 60ml	100/bx



Our open ended catheter tip results in complete nutrition delivery.



Our competitor's rounded closed catheter tip results in pooling of nutrition.

06-0810