



NeoMed® & Medfusion® Enteral Ready Pump Solution Now Featuring the 100 mL Syringe in Legacy and ENFit® Tips

From the syringe pump to the feeding tube, NeoMed delivers an Enteral Safety Solution designed to comply with clinical best practice recommendations.

- Single solution for higher volume nutrition needs of NICU/PICU patients
- Orange or purple faceplate options allow for color coordination with NeoMed's enteral safety system
- Visual differentiation from common parenteral or IV lines
- NeoMed enteral library, now including the newly-characterized 100 mL GravityPro syringe
- Protection from possible off label use from non-characterized syringes



- Fat loss is greater when human milk is administered with a vertical feeding bag, compared to a horizontal syringe*
- Significantly lower priming volumes compared to giving sets used with bags

100 mL NeoConnect Syringe now characterized

- Lower priming volume
- Cost-effective
- Continuity of care for all volume feeds
- Standardization of syringe pump technology
- Offset tip delivers lipids first without tilting pump and voiding the warranty



The Joint Commission

"Use distinctly different pumps for IV applications (rather than using similar pumps for intrathecal and/or epidural applications) to reduce the possibility that an intrathecal medication will accidentally be delivered intravenously and vice versa."

The Joint Commission. "Managing Risk during Transition to New ISO Tubing Connector Standards." Sentinel Event Alert 53 (2014).

* Rayyan, Maissa, et al. "The Fate of Fat: Pre-Exposure Fat Losses during Nasogastric Tube Feeding in Preterm Newborns." *Nutrients*, vol. 7, no. 8, 2015, pp. 6213-6223. doi:10.3390/nu7085279.



Medfusion V6 3500 Enteral Ready Pump with NeoMed Enteral Library

Item #	Product Code
65001	EP-M3500 PUMP

NeoMed Enteral Syringe Software Installation Services	
Product Code	Description
IS-SFW	Cloning Pump
IS-SFW100ML	Cloning Pump; 100 mL Only
IS-SFWSUB	Per Pump for Subsequent Units



Faceplates

Item #	Product Code	Color
65020	EP-MFP ORANGE	Orange
65021	EP-MFP PURPLE	Purple



Pole Clamp

Item #	Product Code
65010	EP-MPC POLE CLAMP




Sterile Clear Syringes with Purple Gradient Markings for Legacy

Item #	Product Code	Size	QTY
11003	PNM-S6EO	6 mL	500
11004	PNM-S12EO	12 mL	500
11005	PNM-S20EO	20 mL	200
11009	PNM-S20EOV	20 mL	200
11006	PNM-S35EO	35 mL	200
11007	PNM-S60EO	60 mL	200
11008	PNM-S100EOV	100 mL	40



Sterile Clear NeoConnect Syringes with Purple Gradient Markings for ENFit

Item #	Product Code	Size	QTY
811003	PNM-S6NC	6 mL	500
811004	PNM-S12NC	12 mL	500
811005	PNM-S20NC	20 mL	200
811009	PNM-S20NCV	20 mL	200
811006	PNM-S35NC	35 mL	200
811007	PNM-S60NC	60 mL	200
811008	PNM-S100NCV	100 mL	40



Sterile Clear Syringes with Orange Gradient Markings for Legacy

Item #	Product Code	Size	QTY
10003	NM-S6EO	6 mL	500
10004	NM-S12EO	12 mL	500
10005	NM-S20EO	20 mL	200
10009	NM-S20EOV	20 mL	200
10006	NM-S35EO	35 mL	200
10007	NM-S60EO	60 mL	200
10008	NM-S100EOV	100 mL	40



Sterile Clear NeoConnect Syringes with Orange Gradient Markings for ENFit

Item #	Product Code	Size	QTY
810003	NM-S6NC	6 mL	500
810004	NM-S12NC	12 mL	500
810005	NM-S20NC	20 mL	200
810009	NM-S20NCV	20 mL	200
810006	NM-S35NC	35 mL	200
810007	NM-S60NC	60 mL	200
810008	NM-S100NCV	100 mL	40



NeoMed offers a variety of purchase and lease options to accommodate most any budgetary need. Please contact your NeoMed Account Manager for additional information.