

October 9, 2018

NEW 100 mL NeoMed® Syringes for use with Medfusion® Syringe Pumps

Dear Customer,

In compliance with the new ISO Standard 80369-3, many syringe manufacturers have launched new enteral syringes, incorporating a new “ENFit®” connection tip.

Smiths Medical has previously validated syringes by model and size for the Medfusion® 3500 and 4000 syringe pumps. Any changes to syringe models and size may have impact on delivery accuracy. The addition of the ENFit® may result in:

- Delivery accuracy confirmation – no modifications were made to the functional design of the syringe, other than adding the new ENFit® connection
- Delivery accuracy impacted – modifications to the syringe with the addition of the new ENFit® connection

New syringe parameters will need to be uploaded on to the pumps with a software update. Once the new syringe software data has been uploaded on the pump the software update will enable all syringes with and without the ENFit® connectors.

NEW 100 mL NeoMed® Syringes

Effective June 1, 2018, the following new NeoMed® syringes are approved for use by Smiths Medical with the Medfusion® 3500 v6 enteral-ready NeoMed® syringe pumps, sold and distributed exclusively by NeoMed®*:

Connector	Syringe Size	Part Number
ENFit® Connector	100 mL	NM-S100NCV
ENFit® Connector	100 mL	PNM-S100NCV
Legacy	100 mL	NM-S100EOV
Legacy	100 mL	PNM-S100EOV

The software update on Medfusion® 3500 v6 pumps includes NeoMed® syringes.

*Note: Characterization data for the 100 mL NeoMed® syringes is approved and available exclusively on Medfusion® 3500 v6 enteral-ready NeoMed® syringe pumps, with the enteral-ready custom software configuration provided exclusively by NeoMed®.

As a reminder, the following NeoMed® syringes were previously approved for use by Smiths Medical with the Medfusion® 3500 pumps (with PharmGuard® v5 or later software) and Medfusion® 4000 pumps (with PharmGuard® v1 or later software).

Syringe Size	NeoMed® 6 mL	NeoMed® 12 mL	NeoMed® 20 mL	NeoMed® 35 mL	NeoMed® 60 mL
NeoMed® ENFit® Syringe Part Numbers	BA-S6NC	BA-S12NC	BA-S20NC	BA-S35NC	BA-S60NC
	BB-S6NC	BB-S12NC	BB-S20NC	BB-S35NC	BB-S60NC
	BC-S6NC	BC-S12NC	BC-S20NC	BC-S35NC	BC-S60NC
	BP-S6NC	BP-S12NC	BP-S20NC	BP-S35NC	BP-S60NC
	NM-S6NC	NM-S12NC	NM-S20NC	NM-S35NC	NM-S60NC
	PNM-S6NC	PNM-S12NC	PNM-S20NC	PNM-S35NC	PNM-S60NC
	BA-S6NC-S	BA-S12NC-S			
		NM-S12NCV	NM-S20NCV		
	PNM-S12NCV	PNM-S20NCV			
NeoMed® Legacy Part Numbers	BA-S6EO	BA-S12EO	BA-S20EO	BA-S35EO	BA-S60EO
	BB-S6EO	BB-S12EO	BB-S20EO	BB-S35EO	BB-S60EO
	BC-S6EO	BC-S12EO	BC-S20EO	BC-S35EO	BC-S60EO
	NM-S6EO	NM-S12EO	NM-S20EO	NM-S35EO	NM-S60EO
	PNM-S6EO	PNM-S12EO	PNM-S20EO	PNM-S35EO	PNM-S60EO
	BA-S6EO-S	BA-S12EO-S			
			NM-S20EOV		
			PNM-S20EOV		

Please contact Smiths Medical U.S. infusion marketing for additional information, or contact your local Smiths Medical infusion account manager at 1.800.258.5361.

Regards,

Hugh L. Springer

Hugh L. Springer
Senior Marketing Manager, Acute Infusion

Tel: +1 763 383 3085 / Mob: +1 651 236 7852
hugh.springer@smiths-medical.com