

Collection Bottles & Transfer Lids

According to the guidelines offered by ASPEN, ADA, HMBANA, and others, proper breast milk collection, transfer, and delivery is essential for adhering to the standard of care. NeoMed's nutrition collection and transfer products were designed to protect and preserve precious nutrition from the home to the hospital setting to meet the specific needs of the neonatal patient.



Robust Gradient Markings

The NeoCollect bottle line offers clearly printed and accurate black gradients in 1 mL increments necessary for precise mixing and compounding.

Concave Floor

Our unique "concave" floor helps pool volume to a narrow area ideal for removal of every single drop of "liquid gold."



Versatile

Universal threads fit most breast pump shields.

ENFit® Compatible Transfer Lid



Oral/Enteral Tip Transfer Lid

Multi-Purpose Transfer Lids

Multi-purpose transfer lids protect contents and facilitate aseptic transfer when filling syringes.

- Fits most standard formula and fortification bottles (Abbott VoluFeed, Mead Johnson, etc.) and are available for both traditional oral or ENFit.
- Conforms to ADA and ASPEN guidelines on aseptic handling of human breast milk.

Orange Enteral/Oral

Description	Code	Item #	QTY
Transfer Lid	TL	90001	200
Collection Bottle 1 oz w/Transfer Lid	SC-1EO	70001	100
Collection Bottle 2 oz w/Transfer Lid	SC-2EO	70002	40

Purple NeoConnect with ENFit Connector

Description	Code	Item #	QTY
Transfer Lid	PTL-NC	891001	200
Collection Bottle 1 oz w/Transfer Lid	PSC-1NC	71001	100
Collection Bottle 2 oz w/Transfer Lid	PSC-2NC	71002	40

NEOBOTTLE®

When delivering nutrition to immunocompromised patients, the NeoBottle* is the only device on the market that offers a closed system barrier to protect HBM from collection to delivery.

The NeoBottle is designed to limit or eliminate breast milk exposure to air, a common limitation with simple open style containers. Ambient air exposure, common with standard screw style lids and bottles, can contribute to the oxidation of milk lipids and an accelerated rate of bacteria colonization.



Leak-Proof Airtight Port

The leak-proof and airtight port may reduce the risk of bacteria colonization which may help to extend the useful life of fresh HBM.

Moveable Floor

NeoMed's engineers created the moveable floor feature to allow the NeoBottle to maintain a closed system barrier as volume is added or removed from the container.



**NeoBottle is only compatible with oral/enteral, not ENFit.*



Amber NeoBottle with Orange Lid

Size	Product Code	Item #	QTY
80 mL	SSNB-AM	12023	60



Amber NeoBottle with Purple Lid

Size	Product Code	Item #	QTY
80 mL	PSSNB-AM	13023	60



Clear NeoBottle with Orange Lid

Size	Product Code	Item #	QTY
80 mL	SSNB-CL	12022	60



Clear NeoBottle with Purple Lid

Size	Product Code	Item #	QTY
80 mL	PSSNB-CL	13022	60